

# HOW TO...

# CONFIGURE NETWORK TASK TRIGGERS?

Last Modified: 06.2018



## **SUMMARY**

D-Lab provides the option to mark tasks or events via remote control. Defined Tasks can be triggered by an event from a simulator or by another third party app via TCP/IP connection.


## **SOFTWARE REQUIREMENTS**

Minimum one D-Lab 3.x Module


## STEP-BY-STEP GUIDE

### Define Tasks and events in D-Lab

1. Define tasks or events in D-Lab in the Task Definitions window.




2. A network name is automatically assigned to each task. This can be changed in D-Lab by double clicking the entry in the Network Name column. Each task must have a unique Network Name (ASCII only).




Name	Network Name	Count	Same Root	Diff Root
Enter Navigation Address	Enter Navigation Address	0	<input type="checkbox"/>	<input type="checkbox"/>
Enter Street	Enter Navigation Address\Enter Street	0		
Enter City	Enter Navigation Address\Enter City	0		
Enter House Number	Enter Navigation Address\Enter House Number	0		
Task Complete	Task Complete	0		

### Add a Network Task Configuration

3. Create a Network Tasks Configuration in the Recording Devices window.  
NOTE: The D-Lab computer and the computer that sends the triggers must be in the same network.



- Establish a TCP connection with the TCP Client (see example code below) and see the feedback in D-Lab



## REQUIREMENTS TO CREATE THE TCP CLIENT

Create a TCP Client and send the network name of the task with the ending "\0" (line 50). It is possible to activate several task at once by sending names in a row (e.g. <task1 network name>\0<task2 network name>\0). For interval tasks, the first network trigger is interpreted as the start and the following one as the end. For single tasks each network trigger marks the task or event.

### Code Example in C#

```
5. using System;
6. using System.Net.Sockets;
7. using System.Text;
8. using System.Threading;

9. namespace DLab3.NetworkTrigger.Demo
10. {
11. ///
12. /// This demo shows how to send Network Task Trigger to D-Lab. If
13. /// running without parameters, it will send the trigger "TEST" to the
14. /// local instance using the default port. When using parameters, host, port,
15. /// trigger name and whether the trigger should be send in a loop can be specified.
16. /// Remember to add and configure the corresponding triggers in D-Lab.
17. ///
18. public class Program
19. {
20. public static void Main(string[] args)
21. {
22. // default values
23. var host = "127.0.0.1";
24. var port = 9000;
25. var trigger = "TEST";
26. var loopDelay = 1000; // value in ms, 0 means no loop

27. // if arguments are provided, use them
28. if (args.Length == 4)
29. {
30. host = args[0];
31. port = int.Parse(args[1]);
32. trigger = args[2];
33. loopDelay = int.Parse(args[3]);
34. }
```

```

34. Console.WriteLine("Demo started. Using these arguments:");
35. Console.WriteLine($" Host: {host}");
36. Console.WriteLine($" Port: {port}");
37. Console.WriteLine($" Trigger: {trigger}");
38. Console.WriteLine($" Loop delay: {loopDelay}");

39. // We need a TCP client. Important: No buffering!
40. using (var tcpClient = new TcpClient { NoDelay = true })
41. {
42. tcpClient.Connect(host, port);
43. Console.WriteLine($"Connected to {host}:{port}.");
44. using (var stream = tcpClient.GetStream())
45. {
46. // If requested we send the trigger every second until a key is pressed
47. do
48. {
49. // We just send the trigger name ASCII encoded with a trailing null terminator
      ("0")
50. var data = Encoding.ASCII.GetBytes(trigger + "\0");
51. stream.Write(data, 0, data.Length);
52. Console.WriteLine($"Sended '{trigger}'.");
53. if (loopDelay > 0) Thread.Sleep(loopDelay);
54. } while (loopDelay > 0 && !Console.KeyAvailable);
55. }
56. }
57. }
58. }
59. }

```

## GET HELP

You can find answers to your questions, as well as instructions and troubleshooting information, in Ergoneers FAQ (<http://www.ergoneers.com/faq>).

In case of other questions or hardware problems please contact [support@ergoneers.com](mailto:support@ergoneers.com).

ERGONEERS GROUP  
WÖHLERWEG 9  
82538 GERETSRIED  
GERMANY

T +49.8171.21624-0  
F +49.8171.21624-11  
INFO@ERGONEERS.COM  
WWW.ERGONEERS.COM

ERGONEERS OF NORTH  
AMERICA, INC.  
111 SW 5TH AVE  
SUITE 3150  
PORTLAND, OR 97204, USA

T +1.503.444.3430  
INFO@ERGONEERS.COM  
WWW.ERGONEERS.COM